



STOP

Do not place full MR64 Clinismart Containers in big yellow bulk bins.



Lock

- 1 Press the red tabs downwards
- 2 Pull the slide latches forward until locked



Label

Always check that the label on top of the MR64 Clinismart container has been completed in full

QC		QA
Hospital name	Hospital ABCD	
Department name	Ward EFGH	
Date opened	01-01-2013	
Date closed	08-01-2013	
Tel: 01388 810310 www.sharpsmart.co.uk		

Place **FULL**
MR64 Clinismart
containers here