



STOP

Do not place full Sharpsmart Containers in big yellow bulk bins.



Lock

- 1 Press the red tabs downwards
- 2 Pull the slide latches forward until locked



Label

Always check that the label on top of the Sharpsmart container has been completed in full

QC		QA
Hospital name	Hospital ABCD	
Department name	Ward EFGH	
Date opened	01-01-2013	
Date closed	08-01-2013	

Tel: 01388 810310 www.sharpsmart.co.uk



Place **FULL**
Sharpsmart
containers here

Sharpsmart
MAKING HEALTHCARE SAFER